

LABXPRT™ v2.0 Label Printer



Part Number	Description
LABXPRT Printer & Accessories	
LABXPRT	v2.0 Printer & Accessories
LABXPRT-KEY	LABXPRT v2.0 LAB LABELING SYSTEM *
LABXPRT-KEY-IDL	LABXPRT v2.0 + IDENTILAB S/WARE-1USER
LABXPRT-KEY-INTL	LABXPRT v2.0 LAB LABELING SYSTEM **
LABXPRT-HC-KEY	HARD CASE FOR LABXPRT V2.0 PRINTER
XPRT-AC (SPARE)	AC ADAPTER for LABXPRT v2.0
XPRT-CABLE-USB	(SPARE) USB PC COMMUNICATIONS CABLE
LABXPRT-SER-USB	SERIAL-TO-USB CABLE ADAPTER (v1 only)

IdentiLab Laboratory Labeling Software

IDL	IDENTILAB LAB ID SOFTWARE-SINGLE USER
IDL-SITE	IDENTILAB LAB ID SOFTWARE-SITE (4+ USERS)





* Ships with: printer, AC power supply, PC cable, XSL-86-461 supply cartridge, CD containing manuals, and driver.

** Ships with: same as LABXPRT-KEY but with international AC power supply.

With Brady's LABXPRT™ v2.0 Label Printer, you have a powerful yet affordable tool at your fingertips! Produce legible, professional looking labels for your lab samples, while drastically reducing the time it actually takes to label. LABXPRT™ printer is now compatible with BradySoft™ Labeling Software.

- Continuous and pre-cut label styles, including round tops, self-laminating labels, and two-label combinations
- Mirror image printing or through-dish viewing of specimen labels
- Recommended for printing 20 to 100 labels at a time
- Continuous and die cut label styles, including tops, self-laminating, and two-label combinations
- Drop, Lock & Go Smart Cartridge combines labels and ribbon for a quick and easy change
- Built-in lab applications, including vial templates, time and date stamp, and auto-serialization for aliquots
- Bar coding capabilities – 2D (when connected to PC) and linear
- 300 dpi print head for sharper resolution on small fonts and bar codes
- Memory to store labels
- 140 Greek and laboratory symbols
- AA battery or AC power operation

LABXPRT™ System Labels



► CENTRIFUGE / EPPENDORF / PCR TUBE LABELS

Diagram	Part Number	Label Width A Inch (mm)	Label Height B Inch (mm)	Write-On Area/ Diameter C Inch (mm)	Labels/ Cartridge	Special Properties
0.2ml PCR TUBE						
⬆	XSL-108-461	0.187 (4.75)	1.500 (38.10)	0.700 (17.78)	500	FNAHSL VS
0.5 TO 0.7ml TUBE						
⬆	XSL-72-461	0.500 (12.70)	1.750 (44.45)	0.750 (19.05)	225	FNAHSL VS
⬆	XSL-123-461	0.375 (9.53)	1.125 (28.58)	0.750 (19.05)	500	FNAHSL VS
⬆	XSL-122-461	0.500 (12.70)	1.125 (28.58)	0.750 (19.05)	500	FNAHSL VS
⊙	XC-375-461	30 ft. (9.1m)	0.375 (9.53)		30 ft. roll	FNAHSL VS
▲	X-130-492	0.825 (20.96)	1.000 (25.40)		275	FNAH L V PG
▲	X-130-499	0.825 (20.96)	0.375 (9.53)		275	FNA V PG
✱	X-118-499	1.000 (25.40)	0.375 (9.53)	0.375 (9.53)	300	FNA V PG
✱	X-120-499	1.000 (25.40)	0.500 (12.70)	0.375 (9.53)	300	FNA V PG
1.0 TO 8.0ml TUBE						
⬆	XSL-74-461	0.500 (12.70)	2.200 (55.88)	0.750 (19.05)	175	FNAHSL VS
⬆	XSL-86-461	0.600 (15.24)	2.625 (66.68)	1.000 (25.40)	150	FNAHSL VS
⬆	XSL-124-461	0.500 (12.70)	1.650 (41.91)	0.750 (19.05)	225	FNAHSL VS
⬆	XSL-126-461	0.600 (15.24)	1.800 (45.72)	1.000 (25.40)	225	FNAHSL VS
⊙	XC-500-461	30 ft. (9.1m)	0.500 (12.70)		30 ft. roll	FNAHSL VS
▲	X-131-492	1.000 (25.40)	0.500 (12.70)		225	FNAH L V PG
▲	X-131-499	1.000 (25.40)	0.500 (12.70)		225	FNA V PG
✱	X-118-499	1.000 (25.40)	0.375 (9.53)	0.375 (9.53)	300	FNA V PG
✱	X-120-499	1.000 (25.40)	0.500 (12.70)	0.375 (9.53)	300	FNA V PG
✱	X-117-499	1.000 (25.40)	0.500 (12.70)	0.440 (11.18)	300	FNA V PG